


Equipment

- X-ray device with an X-ray detector
- Multichannel analyser
- several samples of metal (Zr, Zr, Ni, Ti, Ag, Fe, Cu, Mo or 2 brass, stainless steel, magnet, unknown alloy)
- computer

Setup



Task 1: Calibration of the X-ray device

We set the measurement parameters as the following:

number of channels: 512

negative pulses

amplification factor: -2.5

measurement time: 180s

voltage: 35 kV

current: 1 mA

We also need to calibrate the x-axis. To do this, we fit a Gaussian curve to the K α -lines of Fe and Mo, of which we know the energies ($E_{K\alpha, Fe} \approx 6.4$ keV, $E_{K\alpha, Mo} \approx 17.48$ keV). We record their channels as $\mu_{K\alpha, Fe} \approx 80$, $\sigma_{K\alpha, Fe} \approx 2$, $\mu_{K\alpha, Mo} \approx 216$ and $\sigma_{K\alpha, Mo} \approx 2$, and can thus calibrate the x-axis.

Task 2: Measuring line energies as Gaussian waves

We make the following measurements:

Table 1: Energies of the Gaussian fits for every line

line	μ [keV]	σ [keV]
Zr_{α}	15,84	0,18
Zr_{β}	17,69	0,16
Zn_{α}	8,72	0,19
Zn_{β}	9,66	0,19
Fe_{α}	6,44	0,18
Fe_{β}	7,04	0,32
Mo_{α}	17,49	0,19
Mo_{β}	19,56	0,16
Ni_{α}	7,52	0,19
Ni_{β}	8,28	0,27
Cu_{α}	8,10	0,19
Cu_{β}	8,97	0,25
Ti_{α}	4,61	0,17
Ti_{β}	—	—
Ag_{α}	21,84	0,19
Ag_{β}	24,49	0,23

There was lovely a bump in the curve, and we could not fit a satisfactory Gaussian

Measurement order: Cu, Ag, Zn, Ti, Zr, Ni, Mo, Fe

Task 3: Compositional analysis of alloys

We measure the spectra of our alloy samples, although we do not analyse them numerically. After doing so, we use the programs built-in values for an element's X-ray lines to analyse their composition:

stainless steel: Fe, Ti, V, Mn

brass: Cu, Zn

magnet: Y, Cu, Zn, Ir, Mn

unknown alloy: Cu, Ni, Hg(?)

